

Anti-mouse Langerin/CD207 antibody [923B7] {Alexa-fluor-647} (STJ16102268)

STJ16102268

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-mouse Langerin/CD207 is suitable for use in FACS research applications.
Applications	FACS
Host/Source	Rat
Reactivity	Human/Rat/Mouse/Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	923B7
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-647
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG2a
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	