

Anti-Mannose receptor-like antibody [118F2.02] {Alexa-fluor-488} (STJ16102278)

STJ16102278

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Mannose receptor-like is suitable for use in FACS, Western Blot, Immunohistochemistry and Immunoprecipitation research applications.
Applications	FACS/WB/IHC/IP
Host/Source	Mouse
Reactivity	Human/Sheep

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	118F2.02
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-488
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG2a
Storage	Store at -20°C for up to 1 year.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence