

Anti-pDC/IPC antibody [120G8.04] {Alexa-fluor-488} (STJ16102282)

STJ16102282

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-pDC/IPC is suitable for use in FACS research applications.
Applications	FACS
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	120G8.04
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-488
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1k
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence