

Anti-MIP-3 Alpha/CCL20 antibody [319F6.07] (STJ16102304)
STJ16102304**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-MIP-3 Alpha/CCL20 is suitable for use in ELISA and Immunohistochemistry research applications.
Applications	ELISA/IHC
Host/Source	Mouse
Reactivity	Human/Macaque

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	319F6.07
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG2b
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen
Sequence