

Anti-Neuropilin/CD304 antibody [211H6.01] {Alexa-fluor-647} (STJ16102310)

STJ16102310

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Neuropilin/CD304 is suitable for use in Immunoprecipitation, Immunohistochemistry and FACS research applications.
Applications	IP/IHC-F/FACS
Host/Source	Mouse
Reactivity	Human/Dog/Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	211H6.01
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-647
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage	Store at -20°C for up to 1 year.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence