

Anti-BPAG1 antibody [H10D-Hy-8B] (STJ16102314)

STJ16102314

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Human monoclonal antibody anti-BPAG1 is suitable for use in Western Blot, Immunofluorescence and Immunohistochemistry research applications.
Applications	WB/IF/IHC-F
Host/Source	Human
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	H10D-Hy-8B
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG4
Storage	Store at -20°C for up to 1 year.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence