

Anti-TLR3/CD283 antibody [619F7.06] {Biotin} (STJ16102326)

STJ16102326

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-TLR3/CD283 is suitable for use in Immunohistochemistry, FACS and Western Blot research applications.
Applications	IHC-P/FACS/WB
Host/Source	Mouse
Reactivity	Human/Dog/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	619F7.06
Concentration	0.5 mg/mL
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgM
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	7098
Gene Symbol	TLR3
Uniprot ID	TLR3_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	