

Anti-TLR3/CD283 antibody [713E4.06] {Alexa-fluor-488} (STJ16102333)

STJ16102333

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-TLR3/CD283 is suitable for use in FACS and Immunofluorescence research applications.
Applications	FACS/IF
Host/Source	Mouse
Reactivity	Human/Rat/Mouse/Dog/Sheep/Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	713E4.06
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-488
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	7098
Gene Symbol	TLR3
Uniprot ID	TLR3_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	