

Anti-TLR3/CD283 antibody [716G10.15] {Alexa-fluor-546} (STJ16102339)

STJ16102339

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-TLR3/CD283 is suitable for use in Immunofluorescence research applications.
Applications	IF
Host/Source	Mouse
Reactivity	Human/Rat/Mouse/Dog/Sheep

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	716G10.15
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-546
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgM
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	7098
Gene Symbol	TLR3
Uniprot ID	TLR3_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	