

Anti-TLR3/CD283 antibody [1213F10] {Alexa-fluor-546} (STJ16102344)

STJ16102344

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-TLR3/CD283 is suitable for use in Flow Cytometry, Western Blot, Immunoprecipitation and Immunofluorescence research applications.
Applications	FC/WB/IP/IF
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	1213F10
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-546
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage	Store at -20°C for up to 1 year.
Instruction	

TARGET INFORMATION

Gene ID	7098
Gene Symbol	TLR3
Uniprot ID	TLR3_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	