

**Anti-TLR3/CD283 antibody [1213F10] (STJ16102346)**

STJ16102346

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-TLR3/CD283 is suitable for use in Flow Cytometry, Western Blot, Immunoprecipitation and Immunofluorescence research applications.
<b>Applications</b>	FC/WB/IP/IF
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	1213F10
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	
<b>Formulation</b>	TRIS, PBS, pH7 and antibiotics
<b>Isotype</b>	IgG1
<b>Storage</b>	Store at -20°C for up to 1 year.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">7098</a>
<b>Gene Symbol</b>	<a href="#">TLR3</a>
<b>Uniprot ID</b>	<a href="#">TLR3_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	