

## Anti-TLR3/CD283 antibody [1205C5F1] {Alexa-fluor-488} (STJ16102351)

STJ16102351

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-TLR3/CD283 is suitable for use in Flow Cytometry, Immunohistochemistry and Immunoprecipitation research applications.
<b>Applications</b>	FC/IHC/IP
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	1205C5F1
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Alexa-fluor-488
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	
<b>Formulation</b>	TRIS, PBS, pH7 and antibiotics
<b>Isotype</b>	IgG1
<b>Storage</b>	Store at -20°C for up to 1 year.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">7098</a>
<b>Gene Symbol</b>	<a href="#">TLR3</a>
<b>Uniprot ID</b>	<a href="#">TLR3_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	