

Anti-TLR8/CD288 antibody [303F1.14] {Alexa-fluor-647} (STJ16102357)

STJ16102357

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-TLR8/CD288 is suitable for use in FACS research applications.
Applications	FACS
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	303F1.14
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-647
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG2a
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	170744
Gene Symbol	TLR8
Uniprot ID	TLR8_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	