

Anti-TLR8/CD288 antibody [303F1.14] {Biotin} (STJ16102358)
STJ16102358**GENERAL INFORMATION**

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-TLR8/CD288 is suitable for use in Western Blot research applications. |
| Applications | WB |
| Host/Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|----------------------------------|
| Clonality | Monoclonal |
| Clone ID | 303F1.14 |
| Concentration | 0.5 mg/mL |
| Conjugation | Biotin |
| Purification | Affinity Purified |
| Dilution Range | |
| Formulation | TRIS, PBS, pH7 and antibiotics |
| Isotype | IgG2a |
| Storage Instruction | Store at -20°C for up to 1 year. |

TARGET INFORMATION

| | |
|---------------------------|----------------------------|
| Gene ID | 170744 |
| Gene Symbol | TLR8 |
| Uniprot ID | TLR8_MOUSE |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |