

## Anti-TLR8/CD288 antibody [112H7.15] {Alexa-fluor-647} (STJ16102367)

STJ16102367

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rat monoclonal antibody anti-TLR8/CD288 is suitable for use in FACS research applications.
<b>Applications</b>	FACS
<b>Host/Source</b>	Rat
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	112H7.15
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Alexa-fluor-647
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	
<b>Formulation</b>	TRIS, PBS, pH7 and antibiotics
<b>Isotype</b>	IgG2a
<b>Storage Instruction</b>	Store at -20°C for up to 1 year.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">170744</a>
<b>Gene Symbol</b>	<a href="#">Tlr8</a>
<b>Uniprot ID</b>	<a href="#">TLR8_MOUSE</a>
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	