

Anti-TLR8/CD288 antibody [112H7.15] {Biotin} (STJ16102368)

STJ16102368

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-TLR8/CD288 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	112H7.15
Concentration	0.5 mg/mL
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG2a
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	170744
Gene Symbol	Tlr8
Uniprot ID	TLR8_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	