

Anti-TLR8/CD288 antibody [112H7.15] (STJ16102369)

STJ16102369

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-TLR8/CD288 is suitable for use in Immunohistochemistry, Western Blot and FACS research applications.
Applications	IHC/WB/FACS
Host/Source	Rat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	112H7.15
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG2a
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	170744
Gene Symbol	Tlr8
Uniprot ID	TLR8_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	