

Anti-TLR10/CD290 antibody [25F9] {Alexa-fluor-488} (STJ16102370)

STJ16102370

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-TLR10/CD290 is suitable for use in FACS and Immunohistochemistry research applications.
Applications	FACS/IHC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	25F9
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-488
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	81793
Gene Symbol	TLR10
Uniprot ID	TLR10_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	