

Anti-TLR10/CD290 antibody [9F4] {Alexa-fluor-546} (STJ16102375)

STJ16102375

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-TLR10/CD290 is suitable for use in Immunohistochemistry research applications. |
| Applications | IHC |
| Host/Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|----------------------------------|
| Clonality | Monoclonal |
| Clone ID | 9F4 |
| Concentration | 0.5 mg/mL |
| Conjugation | Alexa-fluor-546 |
| Purification | Affinity Purified |
| Dilution Range | |
| Formulation | TRIS, PBS, pH7 and antibiotics |
| Isotype | IgG1 |
| Storage Instruction | Store at -20°C for up to 1 year. |

TARGET INFORMATION

| | |
|---------------------------|-----------------------------|
| Gene ID | 81793 |
| Gene Symbol | TLR10 |
| Uniprot ID | TLR10_HUMAN |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |