

## Anti-TLR7/CD287 antibody [85D6] {Alexa-fluor-647} (STJ16102385)

STJ16102385

### GENERAL INFORMATION

|                          |                                                                                                                               |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                                            |
| <b>Short Description</b> | Mouse monoclonal antibody anti-TLR7/CD287 is suitable for use in Flow Cytometry and Immunofluorescence research applications. |
| <b>Applications</b>      | FC/IF                                                                                                                         |
| <b>Host/Source</b>       | Mouse                                                                                                                         |
| <b>Reactivity</b>        | Human                                                                                                                         |

### PRODUCT PROPERTIES

|                            |                                  |
|----------------------------|----------------------------------|
| <b>Clonality</b>           | Monoclonal                       |
| <b>Clone ID</b>            | 85D6                             |
| <b>Concentration</b>       | 0.5 mg/mL                        |
| <b>Conjugation</b>         | Alexa-fluor-647                  |
| <b>Purification</b>        | Affinity Purified                |
| <b>Dilution Range</b>      |                                  |
| <b>Formulation</b>         | TRIS, PBS, pH7 and antibiotics   |
| <b>Isotype</b>             | IgG2a                            |
| <b>Storage Instruction</b> | Store at -20°C for up to 1 year. |

### TARGET INFORMATION

|                         |                            |
|-------------------------|----------------------------|
| <b>Gene ID</b>          | <a href="#">170743</a>     |
| <b>Gene Symbol</b>      | <a href="#">Tlr7</a>       |
| <b>Uniprot ID</b>       | <a href="#">TLR7_MOUSE</a> |
| <b>Immunogen</b>        |                            |
| <b>Immunogen Region</b> |                            |
| <b>Specificity</b>      |                            |
| <b>Immunogen</b>        |                            |
| <b>Sequence</b>         |                            |