

Anti-IL1 Alpha antibody [HuMab X3] {Alexa-fluor-546} (STJ16102400)

STJ16102400

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Human monoclonal antibody anti-IL1 Alpha is suitable for use in Universal research applications.
Applications	Neut
Host/Source	Human
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HuMab X3
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-546
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG4/Kappa
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	