

## Anti-Alpha-Synuclein antibody [AS11] {Alexa-fluor-488} (STJ16102409)

STJ16102409

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Alpha-Synuclein is suitable for use in Western Blot, Immunofluorescence and Immunohistochemistry research applications.
<b>Applications</b>	WB/IF/IHC-F
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	AS11
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Alexa-fluor-488
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	
<b>Formulation</b>	TRIS, PBS, pH7 and antibiotics
<b>Isotype</b>	IgG
<b>Storage</b>	Store at -20°C for up to 1 year.
<b>Instruction</b>	

### TARGET INFORMATION

**Gene ID**  
**Gene Symbol**  
**Uniprot ID**  
**Immunogen**  
**Immunogen**  
**Region**  
**Specificity**  
**Immunogen**  
**Sequence**