

## Anti-Alpha-Synuclein antibody [AS11] {Alexa-fluor-647} (STJ16102411)

STJ16102411

### GENERAL INFORMATION

|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Primary antibodies   |
| <b>Short Description</b> | Mouse monoclonal antibody anti-Alpha-Synuclein is suitable for use in Western Blot, Immunofluorescence and Immunohistochemistry research applications. |
| <b>Applications</b>      | WB/IF/IHC-F  |
| <b>Host/Source</b>       | Mouse  |
| <b>Reactivity</b>        | Human  |

### PRODUCT PROPERTIES

|                       |                                  |
|-----------------------|----------------------------------|
| <b>Clonality</b>      | Monoclonal                       |
| <b>Clone ID</b>       | AS11                             |
| <b>Concentration</b>  | 0.5 mg/mL                        |
| <b>Conjugation</b>    | Alexa-fluor-647                  |
| <b>Purification</b>   | Affinity Purified                |
| <b>Dilution Range</b> |                                  |
| <b>Formulation</b>    | TRIS, PBS, pH7 and antibiotics   |
| <b>Isotype</b>        | IgG                              |
| <b>Storage</b>        | Store at -20°C for up to 1 year. |
| <b>Instruction</b>    |                                  |

### TARGET INFORMATION

**Gene ID**  
**Gene Symbol**  
**Uniprot ID**  
**Immunogen**  
**Immunogen**  
**Region**  
**Specificity**  
**Immunogen**  
**Sequence**