

Anti-HIV-1 gp120/160 antibody [113B4.07] {Biotin} (STJ16102414)

STJ16102414

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-HIV-1 gp120/160 is suitable for use in Immunofluorescence, ELISA and FACS research applications.
Applications	IF/ELISA/FACS
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	113B4.07
Concentration	0.5 mg/mL
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen
Sequence