

Anti-HIV-1 gp120/160 antibody [109A2.01] (STJ16102420)
STJ16102420**GENERAL INFORMATION**

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-HIV-1 gp120/160 is suitable for use in FACS and Immunofluorescence research applications. |
| Applications | FACS/IF |
| Host/Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|----------------------------------|
| Clonality | Monoclonal |
| Clone ID | 109A2.01 |
| Concentration | 0.5 mg/mL |
| Conjugation | Unconjugated |
| Purification | Affinity Purified |
| Dilution Range | |
| Formulation | TRIS, PBS, pH7 and antibiotics |
| Isotype | IgG1 |
| Storage Instruction | Store at -20°C for up to 1 year. |

TARGET INFORMATION

| |
|-------------------------|
| Gene ID |
| Gene Symbol |
| Uniprot ID |
| Immunogen |
| Immunogen Region |
| Specificity |
| Immunogen |
| Sequence |