

Anti-HIV-1 gp41 antibody [417D3.08/clone 2] {Alexa-fluor-488} (STJ16102430)
STJ16102430**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-HIV-1 gp41 is suitable for use in ELISA, FACS, Immunofluorescence and Universal research applications.
Applications	ELISA/FACS/IF/Neut
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	417D3.08/clone 2
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-488
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen
Sequence