

## Anti-HIV-1 gp41 antibody [425D11.12/clone 3] {Alexa-fluor-546} (STJ16102431)

STJ16102431

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-HIV-1 gp41 is suitable for use in ELISA, FACS, Immunofluorescence and Universal research applications.
<b>Applications</b>	ELISA/FACS/IF/Neut
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	425D11.12/clone 3
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Alexa-fluor-546
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	
<b>Formulation</b>	TRIS, PBS, pH7 and antibiotics
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	Store at -20°C for up to 1 year.

### TARGET INFORMATION

<b>Gene ID</b>
<b>Gene Symbol</b>
<b>Uniprot ID</b>
<b>Immunogen</b>
<b>Immunogen Region</b>
<b>Specificity</b>
<b>Immunogen</b>
<b>Sequence</b>