

Anti-HIV-1 gp41 antibody [425D11.12/clone 3] {Alexa-fluor-647} (STJ16102433)

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-HIV-1 gp41 is suitable for use in ELISA, FACS, Immunofluorescence and Universal research applications.
Applications	ELISA/FACS/IF/Neut
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	425D11.12/clone 3
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-647
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen
Sequence