

Anti-Myd88 antibody [603E10.06] {Alexa-fluor-488} (STJ16102440)

STJ16102440

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Myd88 is suitable for use in Immunofluorescence research applications.
Applications	IF
Host/Source	Mouse
Reactivity	Human/Dog/Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	603E10.06
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-488
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	17874
Gene Symbol	Myd88
Uniprot ID	MYD88_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	