

Anti-Alpha 4 Beta 7 antibody [204D5.06] {Alexa-fluor-546} (STJ16102445)

STJ16102445

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Alpha 4 Beta 7 is suitable for use in FACS and Immunofluorescence research applications.
Applications	FACS/IF
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	204D5.06
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-546
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	