

Anti-GPIa-IIa antibody [B4F12] {Alexa-fluor-647} (STJ16102473)

STJ16102473

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Human monoclonal antibody anti-GPIa-IIa is suitable for use in FACS and Immunofluorescence research applications.
Applications	FACS/IF
Host/Source	Human
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	B4F12
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-647
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	