

Anti-IL-17A antibody [HAA10E10] {Alexa-fluor-488} (STJ16102482)

STJ16102482

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Human monoclonal antibody anti-IL-17A is suitable for use in Immunofluorescence research applications. |
| Applications | IF |
| Host/Source | Human |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|----------------------------------|
| Clonality | Monoclonal |
| Clone ID | HAA10E10 |
| Concentration | 0.5 mg/mL |
| Conjugation | Alexa-fluor-488 |
| Purification | Affinity Purified |
| Dilution Range | |
| Formulation | TRIS, PBS, pH7 and antibiotics |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for up to 1 year. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | |
| Gene Symbol | |
| Uniprot ID | |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |