

Anti-Chikungunya antibody [E26D9] (STJ16102489)

STJ16102489

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Human monoclonal antibody anti-Chikungunya is suitable for use in Immunofluorescence and Universal research applications.
Applications	IF/Neut
Host/Source	Human
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	E26D9
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1k
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081