

Anti-Chikungunya antibody [E04C1.02] (STJ16102492)

STJ16102492

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Human monoclonal antibody anti-Chikungunya is suitable for use in Immunofluorescence and Universal research applications.
Applications	IF/Neut
Host/Source	Human
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	E04C1.02
Concentration	
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen
Sequence