

Anti-Human IgG1 antibody {Alexa-fluor-546} (STJ16102495)

STJ16102495

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Human polyclonal antibody anti-Human IgG1 is suitable for use in Flow Cytometry and Immunofluorescence research applications.
Applications	FC/IF
Host/Source	Human
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-546
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen
Sequence

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081