

Anti-Human IgG2 antibody {Alexa-fluor-488} (STJ16102498)

STJ16102498

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Human polyclonal antibody anti-Human IgG2 is suitable for use in Flow Cytometry and Immunofluorescence research applications.
Applications	FC/IF
Host/Source	Human
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-488
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG2
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen
Sequence