

Anti-Human IgG3 antibody {Alexa-fluor-647} (STJ16102504)

STJ16102504

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Human polyclonal antibody anti-Human IgG3 is suitable for use in Flow Cytometry and Immunofluorescence research applications.
Applications	FC/IF
Host/Source	Human
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-647
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG3
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen
Sequence