

Anti-Human IgG3 antibody (STJ16102505)

STJ16102505

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Human polyclonal antibody anti-Human IgG3 is suitable for use in Flow Cytometry, Immunofluorescence and Immunohistochemistry research applications.
Applications	FC/IF/IHC
Host/Source	Human
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG3
Storage	Store at -20°C for up to 1 year.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence