

Anti-Human IgG4 antibody {Alexa-fluor-488} (STJ16102506)

STJ16102506

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Human polyclonal antibody anti-Human IgG4 is suitable for use in Flow Cytometry and Immunofluorescence research applications. |
| Applications | FC/IF |
| Host/Source | Human |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|----------------------------------|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.5 mg/mL |
| Conjugation | Alexa-fluor-488 |
| Purification | Affinity Purified |
| Dilution Range | |
| Formulation | TRIS, PBS, pH7 and antibiotics |
| Isotype | IgG4 |
| Storage Instruction | Store at -20°C for up to 1 year. |

TARGET INFORMATION

| |
|-------------------------|
| Gene ID |
| Gene Symbol |
| Uniprot ID |
| Immunogen |
| Immunogen Region |
| Specificity |
| Immunogen |
| Sequence |