

Anti-Human IgG4 antibody (STJ16102509)

STJ16102509

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Human polyclonal antibody anti-Human IgG4 is suitable for use in Flow Cytometry, Immunofluorescence and Immunohistochemistry research applications. |
| Applications | FC/IF/IHC |
| Host/Source | Human |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|-----------------------|----------------------------------|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.5 mg/mL |
| Conjugation | Unconjugated |
| Purification | Affinity Purified |
| Dilution Range | |
| Formulation | TRIS, PBS, pH7 and antibiotics |
| Isotype | IgG4 |
| Storage | Store at -20°C for up to 1 year. |
| Instruction | |

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence