

Anti-Mouse IgG2a antibody (STJ16102526)

STJ16102526

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse polyclonal antibody anti-Mouse IgG2a is suitable for use in FACS, Immunofluorescence, Immunohistochemistry and Immunohistochemistry research applications.
Applications	FACS/IF/IHC-F/IHC-P
Host/Source	Mouse
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG2a
Storage	Store at -20°C for up to 1 year.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence