

Anti-Rat IgG2a antibody {Alexa-fluor-647} (STJ16102545)
STJ16102545**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Rat polyclonal antibody anti-Rat IgG2a is suitable for use in FACS and Immunofluorescence research applications.
Applications	FACS/IF
Host/Source	Rat
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-647
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG2a
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	