

Anti-Rat IgG2b antibody (STJ16102546)

STJ16102546

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat polyclonal antibody anti-Rat IgG2b is suitable for use in FACS, Immunofluorescence, Immunohistochemistry and Immunohistochemistry research applications.
Applications	FACS/IF/IHC-F/IHC-P
Host/Source	Rat
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG2b
Storage	Store at -20°C for up to 1 year.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence