

Anti-Rat IgG2b antibody {Alexa-fluor-546} (STJ16102548)
STJ16102548**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Rat polyclonal antibody anti-Rat IgG2b is suitable for use in FACS and Immunofluorescence research applications.
Applications	FACS/IF
Host/Source	Rat
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-546
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG2b
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence