

**Anti-Rat IgG3 antibody (STJ16102550)**

STJ16102550

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rat polyclonal antibody anti-Rat IgG3 is suitable for use in FACS, Immunofluorescence, Immunohistochemistry and Immunohistochemistry research applications.
<b>Applications</b>	FACS/IF/IHC-F/IHC-P
<b>Host/Source</b>	Rat
<b>Reactivity</b>	Rat

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	
<b>Formulation</b>	TRIS, PBS, pH7 and antibiotics
<b>Isotype</b>	IgG3
<b>Storage</b>	Store at -20°C for up to 1 year.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>
<b>Gene Symbol</b>
<b>Uniprot ID</b>
<b>Immunogen</b>
<b>Immunogen</b>
<b>Region</b>
<b>Specificity</b>
<b>Immunogen</b>
<b>Sequence</b>