

Anti-Mouse IgG3 antibody {Alexa-fluor-647} (STJ16102553)

STJ16102553

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat polyclonal antibody anti-Mouse IgG3 is suitable for use in FACS and Immunofluorescence research applications.
Applications	FACS/IF
Host/Source	Rat
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-647
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG3
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence