

Anti-Estrogen receptor antibody [AER-320] (STJ16102572)

STJ16102572

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Estrogen receptor is suitable for use in Immunohistochemistry research applications.
Applications	IHC
Host/Source	Mouse
Reactivity	Cow/Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	AER-320
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	Immunohistology (1-2 µg/ml for 30 min at RT; an appropriate antigen retrieval method for staining of formalin-fixed tissues has not been established to date).
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG2ak
Storage	Store for up to 1 year at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	2852
Gene Symbol	GPB1
Uniprot ID	GPB1_HUMAN
Immunogen	A BALB/c mouse was immunized with estrogen receptor from calf uterus.Fusion partner: Sp2/0.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	