

## Anti-Catenin p120, delta-1-TNND1, pTyr96 antibody [EBS-CA-011] (STJ16102583)

STJ16102583

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Catenin p120, delta-1-TNND1, pTyr96 is suitable for use in Western Blot and Immunohistochemistry research applications.
<b>Applications</b>	WB/IHC
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human/Mouse/Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	EBS-CA-011
<b>Concentration</b>	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purified from tissue culture.
<b>Dilution Range</b>	Immunoblotting (1-2 µg/ml). Immunohistology (1-2 µg/ml for 30 min at RT; an appropriate antigen retrieval method for staining of formalin-fixed tissues has not been established to date).
<b>Formulation</b>	PBS with 0.02% Sodium Azide.
<b>Isotype</b>	IgG1k
<b>Storage</b>	Store for up to 1 year at 2-8°C upon receipt.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	A BALB/c mouse was immunized with Mouse p120 Catenin (pY96).Fusion partner: SP2/0.
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	