

Anti-CD36 antibody [EBS-CD-024] (STJ16102587)

STJ16102587

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD36 is suitable for use in Flow Cytometry, Immunofluorescence and Immunohistochemistry research applications.
Applications	FC/IF/IHC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EBS-CD-024
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	Flow cytometry (1-2 µg/million cells in 0, 1 ml for 30 min, at 4°C). Immunofluorescence (1-2 µg/ml). Immunohistology (1-2 µg/ml for 30 min at RT; no antigen retrieval procedure is necessary for staining of formalin-fixed tissues).
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG1k
Storage	Store for up to 1 year at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	948
Gene Symbol	CD36
Uniprot ID	CD36_HUMAN
Immunogen	A BALB/c mouse was immunized with gp88 purified from human platelets. Fusion partner: SP2/0.
Immunogen Region	
Specificity	
Immunogen Sequence	