

Anti-CD176/TF Antigen antibody [EBS-CD-054] (STJ16102593)

STJ16102593

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD176/TF antigen is suitable for use in ELISA, Western Blot and Immunohistochemistry research applications.
Applications	ELISA/WB/IHC
Host/Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EBS-CD-054
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	ELISA (solid phase: not known; tracer: 0, 001-100 µg/ml for 30 min at RT). Immunoblotting (1-2 µg/ml). Immunohistology (1-2 µg/ml for 30 min at RT).
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgMκ
Storage	Store for up to 1 year at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	A BALB/c mouse was immunized with neuraminidase-treated human red blood cells. Fusion partner: X63-AG8.653.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	